

By J.R. Garbarino

Arsenic Field-Speciation Method Checklist and Worksheet (Schedule 1729)

Site ID: _____ Date/Time: _____
 Site Name: _____ Collector: _____

- _ Filter sample using 0.45- μ m capsule filter; do not clean media with acid
- _ SPE Cartridge Conditioning: 2 mL methanol and 10 mL deionized water – pH _____

or

- _ SPE Cartridge Conditioning: 5 mL 1.7 M acetic acid and 10 mL deionized water – pH _____

- _ Estimated anion milliequivalents in 10 ml sample: _____ meq

$$C_a = [\text{HCO}_3^- \times 1.6(10^{-4})] + [\text{Cl}^- \times 2.8(10^{-4})] + [\text{NO}_3^- \times 7.1(10^{-4})] + [\text{HPO}_4^{-2} \times 2.1(10^{-4})] + [\text{SO}_4^{-2} \times 2.1(10^{-4})]$$

where: C_a = anion concentration in milliequivalents in 10-mL sample

HCO_3^- = bicarbonate concentration in mg/L as HCO_3^-

Cl^- = chloride concentration in mg/L as Cl^-

NO_3^- = nitrate concentration in mg/L as N

HPO_4^{-2} = phosphate concentration in mg/L as HPO_4^{-2}

SO_4^{-2} = sulfate concentration in mg/L as SO_4^{-2}

- _ Diluent (blank water) volume added to tube A to give <0.1 anion milliequivalents in 10-mL sample _____ mL
- _ EDTA preservative added to tube A _____ μ L

(100 μ L or the volume calculated below, whichever is greater)

(Note: if sample is diluted, account for dilution factor when using the following equation):

$$V_{\text{EDTA}} = 4.0(10^6) \times ([\text{Al} \times 3.7(10^{-10})] + [\text{Fe} \times 1.8(10^{-10})] + [\text{Mn} \times 1.8(10^{-10})] + [\text{Ca} \times 2.5(10^{-7})] + [\text{Mg} \times 4.1(10^{-7})] + [\text{Sr} \times 1.1(10^{-10})])$$

where: V_{EDTA} = microliters of 250-mM EDTA required per sample

Al = dissolved aluminum concentration, in μ g/L as Al

Fe = dissolved iron concentration, in μ g/L as Fe

Mn = dissolved manganese concentration, in μ g/L as Mn

Ca = dissolved Ca concentration, in mg/L as Ca

Mg = dissolved magnesium concentration, in mg/L as Mg

Sr = dissolved strontium concentration, in μ g/L as Sr

- _ Sample volume added to tube A (10.0 \pm 0.2 or 10.0 mL minus diluent volume) _____ mL
- _ Spike Solution lot number _____ F
 Solution concentration _____ μ g/L
 Volume added _____ μ L
- _ Push all the sample through cartridge at 1 to 2 drops per second.
- _ Collect eluate in tube B; cap tightly.
- _ Write Station ID, date, and time on the cartridge, tube B, and shipping container.
- _ Place tube B, cartridge (inside its shipping container), and copy of worksheet inside the kit bag. Label bag with station ID, date, and time.
- _ Maintain at 4 $^{\circ}$ C. Ship chilled sample and a copy of the worksheet to the NWQL within 14 days of collection.

COMMENTS:

Figure 5-5. Worksheet for field-speciation method to determine arsenic species in water samples.

Arsenic Laboratory-Speciation Methods Checklist and Worksheet
(Schedule 1730, 1731, or 1732)

Site ID: _____ Date/Time: _____
Site Name: _____ Collector: _____
Lab schedule requested: _____

Laboratory-speciation methods Schedules 1730, 1731, or 1732

- Filter sample using 0.45- μ m disposable capsule filter;
do not clean media with acid
- EDTA preservative added to opaque bottle _____ μ L
(100 μ L or the volume calculated below, whichever is greater)

$$V_{\text{EDTA}} = 4.0(10^{-6}) \times ([\text{Al} \times 3.7(10^{-10})] + [\text{Fe} \times 1.8(10^{-10})] + [\text{Mn} \times 1.8(10^{-10})] + [\text{Ca} \times 2.5(10^{-7})] + [\text{Mg} \times 4.1(10^{-7})] + [\text{Sr} \times 1.1(10^{-10})])$$

where: V_{EDTA} = microliters of 250-mM EDTA required per sample
Al = dissolved aluminum concentration, in $\mu\text{g/L}$ as Al
Fe = dissolved iron concentration, in $\mu\text{g/L}$ as Fe
Mn = dissolved manganese concentration, in $\mu\text{g/L}$ as Mn
Ca = dissolved Ca concentration, in mg/L as Ca
Mg = dissolved magnesium concentration, in mg/L as Mg
Sr = dissolved strontium concentration, in $\mu\text{g/L}$ as Sr

- Spike Solution lot number _____ **L or F**
(circle one)
Solution concentration _____ $\mu\text{g/L}$
Volume added _____ μL

Spiked or unspiked sample volume completely fill bottle, but do not overflow.

Note "8-mL" opaque bottle contains 11.5 ± 0.1 mL when completely full: 11.5 mL

- Write Station ID, date, and time on bottle.
- Maintain at 4 °C. Ship chilled sample and a copy of the worksheet to the NWQL within 14 days of collection.

COMMENTS:

Figure 5-6. Worksheet for laboratory-speciation methods to determine arsenic species in water samples.