

Report for 2004MO31B: Use of Excitation/Emission Matrix Fluorescence

There are no reported publications resulting from this project.

Report Follows

An annual report submitted to Missouri Water Resources Research Center/USGS

Use of Excitation/Emission Matrix Fluorescence for Wastewater Sources Identification

Dr. Baolin Deng
Department of Civil and Environmental Engineering
University of Missouri-Columbia
Columbia, MO 65211

ABSTRACT

Water contamination can come from different sources, including wastewater treatment effluent, landfill leachate, agricultural runoff, and a wide variety of industrial sources. Identifying and differentiating these sources are needed to effectively manage the water pollution issues and enforce the environmental regulation. The objective of this project is to establish a fluorescence excitation/emission matrix database for different types of water and to develop a mathematical procedure/computer model to identify the source of unknown water samples. The research tasks include: (i) collections of substantial numbers of water samples from various known sources, (ii) analysis of water samples for excitation/emission spectra, and (iii) development of a mathematical/computer model that can be used to identify different sources of water. In this project period, we have made significant progress on tasks (i) and (ii), including procurement, installation, and testing of the a Hitachi F-4500 fluorescence spectrophotometer; analysis of pure water samples and individual chemical dissolved in water; and analysis of 60 Missouri water samples from various lakes and water supplies.

Introduction

Molecular fluorescence spectroscopy is an analytical technique based on emission of fluorescence from a molecule upon excitation by irradiation. It is a rapid and inexpensive way of wastewater identification and quantification, in particular, when excitation/emission matrix spectra (EEMS) in a wide range of wavelength are collected. Compounds containing aromatic rings normally give the most intense molecular fluorescence emission, including benzene, toluene, phenol, quinoline, and many fused ring aromatic compounds. A small number of aliphatic and alicyclic carbonyl compounds with highly conjugated double-bond structures also fluorescence. Since fluorescence spectroscopy is much more sensitive than absorption spectroscopy, it is particularly suitable for trace amounts of contaminant detection.

The objective of this project is to establish a fluorescence EEMS database for different types of water in the state of Missouri and develop a mathematical procedure/computer model to identify the source of unknown water samples. Different types of waters have their own fluorescence features as they contain various fluorescence species. The fluorescence characteristics of river water (Yan et al., 2000), seawater (Coble et al., 1990; Coble, 1996), groundwater (Baker and Genry, 1999; Baker and Lamont-Black, 2001), sewage, and sewage-impacted river water (Reynolds and Ahmad, 1997; Baker, 2001, 2002) have been documented separately. However, by now, there is no database that provides comprehensive fluorescence spectra for various types of water. This research aims to fill the data gap by building a much more comprehensive database critical for the source identification of the unknown wastewater samples.

Experimental Section

1. Instrument Description and Parameters for Data Collection

We have acquired a Hitachi F-4500 fluorescence spectrophotometer for the current study. This is a research grade, easy to use fluorescence spectrophotometer with its FL Solutions[®] software for rapid data acquisition (2400nm/min-30,000 nm/min).

Instrument Features (Courtesy of Hitachi)

- Three-dimensional contour plots provide maximum spectral information and can be acquired and displayed in less than one minute.
- Phosphorescence is built-in as standard and allows lifetime measurements from 1 ms to address the widest range of samples.
- In addition to fluorescence and phosphorescence measurements, luminescence is included as standard. The high energy throughput and signal-to noise of the system allow measurements of chemiluminescent and bioluminescent compounds.
- The intense 150 watt xenon source provides maximum light energy over the entire 200-900 nm wavelength range of the spectrophotometer. A rhodamine B quantum counter is used to determine the correction factors for true instrument response to assure spectral accuracy.

- Unique horizontal beam geometry of the excitation and emission beams increases sensitivity and requires only 0.6 mL of sample in a standard 10mm cuvette.
- Automatic shutter control protects photo sensitive samples from degradation.

The following parameters, as illustrated in Table 1, were set before using the instrument for measurement of excitation and emission spectra. The excitation wavelength was set from 200 to 400 nm and emission wavelength from 250 to 550 nm on 12000nm/min and 30000nm/min scan speed.

Table 1: Calibration set for scanning

| Information about sample/instrument | | Calibration set for instrument | |
|-------------------------------------|-------------------------------|--------------------------------|--------------|
| Report Date: | 5/19/2005 22:11 | EX End WL: | 400.0 nm |
| Sample: | Watkins Mill(#74) | EX Sampling Interval: | 2.0 nm |
| File name: | Missouri Water Samples004.FD3 | EM Start WL: | 250.0 nm |
| Run Date: | 3/11/2005 9:36 | EM End WL: | 550.0 nm |
| Operator: | Owner | EM Sampling Interval: | 3.0 nm |
| Comment: | Collected @ 08/09/04 | Scan speed: | 12000 nm/min |
| Instrument Model: | F-4500 FL Spectrophotometer | EX Slit: | 5.0 nm |
| Serial Number: | 1606-005 | EM Slit: | 10.0 nm |
| ROM Version: | 4000 05 | PMT Voltage: | 700 V |
| Measurement type: | 3-D Scan | Response: | Auto |
| Data mode: | Fluorescence | Corrected spectra: | On |
| EX Start WL: | 200.0 nm | Shutter control: | On |

2. Analysis

To establish a database for water source identification, we need to collect water samples from various known sources and analyze the samples for excitation/emission matrix spectra. Up to now, 60 water samples collected from major lakes and water treatment plants and were analyzed. The Milli-Q deionized water and tap water, as well as solutions of individual chemicals with various functional groups were also scanned for excitation and emission spectra to understand the fluorescence characteristics and pattern and intensities. Fluorescence data were collected in the excitation wavelength from 200 to 400 nm and emission wavelength from 250 to 550 nm with 2 and 3 nm intervals respectively.

1. Tap & D.I Water Samples:

Several deionized water samples and tap water from University of Missouri Campus were scanned for excitation and emission matrix spectra. Following are the 3-D scan of each of them. The former shows no fluorescence and latter has some degree of intensity and fluorescence possible because of natural organic matter.




Figure 1




Figure 2

2. Water Samples for Individual Chemicals:

Eighteen chemicals, as listed in Table 2, were scanned at a concentration of 10 ppm for each.

Table 2: List of chemicals and water samples

| S.N | Sample Name | Volume (ml) | Conc. (g/l) | Weight (g) |
|-----|---|-------------|-------------|------------|
| 1 | Benzoic Acid | 50 | 5 | 0.25 |
| 2 | Glutamic Acid | 50 | 5 | 0.25 |
| 3 | Hydroxybenzoic Acid | 50 | 5 | 0.25 |
| 4 | Phenol | 50 | 5 | 0.25 |
| 5 | D-Glucose | 50 | 5 | 0.25 |
| 6 | 4-amino- Pyrine | 50 | 5 | 0.25 |
| 7 | 1,2 cyclohexylene-dinitrilo tetracetic acid | 50 | 5 | 0.25 |
| 8 | Cosmic acid | 50 | 5 | 0.25 |
| 9 | 5'-5''-Dibromo-O- creasolesulfonephthalein | 50 | 5 | 0.25 |
| 10 | Dimethylamino azobenzene | 50 | 5 | 0.25 |
| 11 | 1-5 Diphenylcarbazide | 50 | 5 | 0.25 |
| 12 | ADA,N(2-Acetamido) Iminodiacetic Acid | 50 | 5 | 0.25 |
| 13 | Oxalic Acid | 50 | 5 | 0.25 |
| 14 | 3-(4-Morpholino)-Propanesulfonic Acid | 50 | 5 | 0.25 |
| 15 | 2-(4-Morpholino)-Ethanesulfonic Acid | 50 | 5 | 0.25 |
| 16 | 1,10 Phenanthroline | 50 | 5 | 0.25 |
| 17 | Peptone | 50 | 5 | 0.25 |
| 18 | 2-(2-Thiazolylazo)-p-Cresol | 50 | 5 | 0.25 |

The data collected could be visualized by 3-D contour map and Bird's eye view, as illustrated by Figure 3-6 for hydroxybenzoic acid and D-glucose.




Figure 3: Contour and spectrum of hydroxybenzoic acid




Figure 4: Bird's eye view of Hydroxybenzoic acid




Figure 5 D-Glucose




Figure 6: Bird's eye view of D-Glucose

3. Missouri Water Samples

A total of sixty (60) samples were collected through Missouri Water Resources Research Center and scanned for excitation and emission spectra. The samples represent the different rivers, lakes, creeks in the State of Missouri. Table 4 lists the samples scanned for the fluorescence intensity and some EEMS are illustrated in Figures 7-14. .

Table 4: A list of water samples collected for fluorescence analysis

| S.N | Collected Date | CE Bottle # | Sample Name | Sample# |
|------------|-----------------------|------------------------|-------------------------|----------------|
| 1 | 7/27/2004 | 7/27-1 | Miller Community Lake | 35 |
| 2 | 7/27/2004 | 7/27-3 | Clearwater Lake | 36 |
| 3 | 7/27/2004 | 7/27-5 | Wapapello Lake | 30 |
| 4 | 7/27/2004 | 7/27-7 | Lake Girardeau | 28 |
| 5 | 7/27/2004 | 7/27-9 | Boutin Lake | 27 |
| 6 | 7/27/2004 | 7/27-11 | Bella Vista Lake | 26 |
| 7 | 7/27/2004 | 7/27-13 | Fredericktown City Lake | 24 |
| 8 | 7/28/2004 | 7/28-15 | Anne | 22 |
| 9 | 7/28/2004 | 7/28-17 | Capri | 18 |
| 10 | 7/28/2004 | 7/28-19 | Shane | 112 |
| 11 | 7/28/2004 | 7/28-21 | Bismark | 186 |
| 12 | 7/28/2004 | 7/28-23 | Council Bluff | 39 |
| 13 | 8/3/2004 | 8/3-1 | Monroe Lake | 60 |
| 14 | 8/3/2004 | 8/3-3 | Lake Hunnewell | 115 |
| 15 | 8/3/2004 | 8/3-5 | Henry Sever | 114 |
| 16 | 8/3/2004 | 8/3-7 | LaBelle City | 129 |
| 17 | 8/3/2004 | 8/3-9 | Deer Ridge | 57 |
| 18 | 8/3/2004 | 8/3-11 | Fox Valley | 182 |
| 19 | 8/3/2004 | 8/3-13 | Lancaster City | 132 |
| 20 | 8/4/2004 | 8/4-15 | LaPlatta Lake | 130 |
| 21 | 8/4/2004 | 8/4-17 | Longbranch | 48 |
| 22 | 8/4/2004 | 8/4-19 | Thomas Hill Reservoir | 46 |
| 23 | 8/4/2004 | 8/4-21 | Tri-CityCommunity Lake | 42 |
| 24 | 8/4/2004 | 8/4-23 | Little Dixie | 117 |
| 25 | 8/9/2004 | 8/9-1 | Pape Lake | 113 |
| 26 | 8/9/2004 | 8/9-2 | Higginsville | 121 |
| 27 | 8/9/2004 | 8/9-3 | Hazel Hill | 183 |
| 28 | 8/9/2004 | 8/9-4 | Odessa Lake | 64 |
| 29 | 8/9/2004 | 8/9-5 | Watkins Mill | 74 |
| 30 | 8/9/2004 | 8/9-6 | Smithville Lake | 72 |
| 31 | 8/9/2004 | 8/9-7 | Belcher Branch | 187 |
| 32 | 8/10/2004 | 8/10-8 | Raintree | 118 |
| 33 | 8/10/2004 | 8/10-9(not filtrated) | Holden City | 157 |
| 34 | 8/10/2004 | 8/10-10(not filtrated) | North Lake | 70 |
| 35 | 8/10/2004 | 8/10-11(not filtrated) | Butler City | 189 |
| 36 | 8/10/2004 | 8/10-12(not filtrated) | Atkinson | 91 |
| 37 | 8/10/2004 | 8/10-13(not filtrated) | Bushwacker | 161 |
| 38 | 8/10/2004 | 8/10-14(not filtrated) | Stockton | 93 |
| 39 | 8/10/2004 | 8/10-15(not filtrated) | McDaniel | 95 |
| 40 | 8/10/2004 | 8/10-16(not filtrated) | Fellows | 96 |
| 41 | 8/11/2004 | 8/11-17(not filtrated) | Pomme De Terre | 92 |
| 42 | 8/11/2004 | 8/11-18(not filtrated) | Truman | 89 |

| | | | | |
|----|-----------|------------------------|--------------------|-----|
| 43 | 8/11/2004 | 8/11-19(not filtrated) | Lake of the Ozarks | 149 |
| 44 | 8/11/2004 | 8/11-20(not filtrated) | Manito | 184 |
| 45 | 8/11/2004 | 8/11-21(not filtrated) | Binder | 13 |
| 46 | 8/17/2004 | 8/17-1 | Sterling | 119 |
| 47 | 8/17/2004 | 8/17-2 | Nehai Tonkyea | 45 |
| 48 | 8/17/2004 | 8/17-3 | Marceline-2 | 88 |
| 49 | 8/17/2004 | 8/17-4 | Marceline-1 | 139 |
| 50 | 8/17/2004 | 8/17-5 | Brookfield | 87 |
| 51 | 8/17/2004 | 8/17-6 | Hazel Creek | 131 |
| 52 | 8/17/2004 | 8/17-7 | Forest | 133 |
| 53 | 8/17/2004 | 8/17-8 | Green City | 137 |
| 54 | 8/18/2004 | 8/18-9 | Marie | 85 |
| 55 | 8/18/2004 | 8/18-10 | Paho | 84 |
| 56 | 8/18/2004 | 8/18-11 | New Bethany | 83 |
| 57 | 8/18/2004 | 8/18-12 | Harrison Co. | 185 |
| 58 | 8/18/2004 | 8/18-13 | Mozingo | 181 |
| 59 | 8/18/2004 | 8/18-14 | Nodaway | 179 |
| 60 | 8/18/2004 | 8/18-15 | Bilby Ranch | 180 |




Figure 7




Figure 8




Figure 9




Figure 10




Figure 11




Figure 12




Figure 13




Figure 14

Summary

We have made significant progress in the first stage of this research, including (i) procurement, installation, and testing of the a Hitachi F-4500 fluorescence spectrophotometer; (ii) analysis of pure water samples and individual chemical dissolved in water; and (iii) analysis of 60 Missouri water samples from various lakes and water supplies. Different chemicals show substantially

different characteristics of fluorescence spectra. Water samples also show different characteristics. Further research is under way to establish the relationship between the fluorescence features of individual chemicals and natural water samples, and between fluorescence features of natural water samples and water quality parameters such as NDMA formation potentials.

Literature cited

- (1) Baker, A. Fluorescence excitation-emission matrix characterization of some sewage-impacted rivers. *Environ. Sci. Technol.*, 2001, 35, 948-953.
- (2) Baker, A. Fluorescence excitation-emission matrix characterization of river waters impacted by tissue mill effluent. *Environ. Sci. Technol.* 2002, 36, 1377-1382.
- (3) Baker, A.; Genry, D. Fluorescence wavelength and intensity variations of cave waters. *Journal of Hydrology*, 1999, 217, 19-34.
- (4) Baker, A.; Lamont-Black, J. Fluorescence of dissolved organic matter as a natural tracer of ground water. *Ground Water*, 2001, 39(5), 745-750.
- (5) Coble, P.G. Characterization of marine and terrestrial DOM in seawater using excitation-emission matrix spectroscopy. *Marine Chemistry*, 1996, 51,325-346.
- (6) Coble, P.G.; Green, S.A.; Blough, N.V.; Gagosian, R.B. Characterization of dissolved organic matter in the Black Sea by fluorescence spectroscopy. *Nature*, 1990, 348, 432-435.
- (7) Reynolds, D.M.; Ahmad, S.R. Rapid and direct determination of wastewater BOD values using a fluorescence technique. *Water Research*, 1997, 31, 2012-2018.
- (8) Yan, Y; Li, H.; Myrick, M.L. Fluorescence fingerprint of waters: excitation-emission matrix spectroscopy as a tracking tool. *Applied spectroscopy*, 2000, 54(10),1539-1542.